


Purified anti-human Podoplanin Antibody



SeqIF™ (sequential immunofluorescence) staining on COMET™ of Purified anti-Podoplanin (clone NC-08, yellow) on formalin-fixed paraffin-embedded human skin tissue at 5 µg/mL. Alexa Fluor™ Plus 555 Goat anti-Rat IgG antibody (Lunaphore, Cat. No. DR555RT) was used as a secondary antibody. Nuclei were counterstained with DAPI (blue). Tissue underwent an all-in-one dewaxing and antigen retrieval preprocessing.