


Recombinant Mouse Sonic Hedgehog (Shh) (Cys25/(Ile-Ile)-Gly198) (carrier-free)



Recombinant Mouse Sonic Hedgehog (SHH) (Cys25/(Ile-Ile)-Gly198) induces C3H10T1/2 cell differentiation in the presence of BMP9 (20 ng/mL) as measured by alkaline phosphatase production. Alkaline phosphatase activity was measured using p-nitrophenyl phosphate. BioLegend's protein was compared side-by-side to a leading competitor's equivalent product.