

LEGENDplex™ Mouse CCL9 (MIP-1γ) Capture Bead B9, 13X

Catalog# / Size	741079 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	CCF18, Macrophage inflammatory protein (MIP)-related protein-2 (MRP-2)
Description	<p>The beads are intended for use with the following Mix and Match assay reagents.</p> <p>Note: for serum/plasma samples, MCP-2 and MIP-1γ cannot be mixed and matched with other targets and need to be measured separately due to high dilution requirement.</p> <p>741070 LEGENDplex™ MU Proinflam. Chemokine Panel 2 Standard 741145 LEGENDplex™ Buffer Set O 741069 LEGENDplex™ MU Proinflam. Chemokine Panel 2 Detection Antibodies 740377 LEGENDplex™ Filter Plate OR 740379 LEGENDplex™ V-bottom Plate</p>

Product Details

Reactivity	Mouse
Application	Multiplex
	Learn more about LEGENDplex™ at biolegend.com/legendplex
	Download the LEGENDplex™ software here .

Antigen Details

Gene ID	20308
----------------	-----------------------

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587