

LEGENDplex™ Carboxyl Beads A5

Catalog# / Size	740168 / 1 vial
Regulatory Status	RUO
Description	LEGENDplex™ Carboxyl Beads are fluorescent particles designed for the multiple assay development using flow cytometers. Beads include two sizes, i.e., 5 µm and 7 µm, designated as A or B beads, and each bead set has different levels of APC fluorescence, allowing them to be clearly distinguished from each other. The functionalized carboxyl group can be used for the conjugation of antibody and streptavidin.

Product Details

Formulation	Each vial contains 10 ⁷ particles (10 million) per mL.
Storage & Handling	Store between 2°C and 8°C.
Application	Multiplex Learn more about LEGENDplex™ at biolegend.com/legendplex Download the LEGENDplex™ software here .
Application Notes	LEGENDplex™ Carboxyl Beads are intended for do-it-yourself conjugations to create your own assays. Please refer to the conjugation protocol here .


Antigen Details

Gene ID	NA
----------------	----

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

Product Data



Beads A: Six sets of smaller beads separated by FL4 (APC) Channel on FACSCalibur™

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally

similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587