

## LEGENDplex™ Mouse CXCL5 (LIX) Capture Bead B7, 13X

<b>Catalog# / Size</b>	740097 / 100 tests (270 µL)
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	CXCL5
<b>Description</b>	LEGENDplex™ Mouse Proinflammatory Chemokine Mix and Match Panel Capture Beads are intended for use with the following reagents:  740074 (LEGENDplex™ Mouse Proinflammatory Chemokine Detection Antibodies) 740370 (LEGENDplex™ Mouse Proinflammatory Chemokine Standard) 740368 (LEGENDplex™ Buffer Set A)

### Product Details

---

<b>Verified Reactivity</b>	Mouse
<b>Application</b>	<a href="#">Multiplex</a>  Learn more about LEGENDplex™ at <a href="http://biolegend.com/legendplex">biolegend.com/legendplex</a>  Download the <a href="#">LEGENDplex™ software here</a> .

<b>Product Citations</b>	1. Popescu M, <i>et al.</i> 2019. J Immunol. 203:2415. <a href="#">PubMed</a>
--------------------------	---

### Antigen Details

---

<b>Biology Area</b>	Immunology, Innate Immunity
<b>Molecular Family</b>	Cytokines/Chemokines
<b>Gene ID</b>	<a href="#">20311</a>

### Related Protocols

---

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

