

LEGENDplex™ Human BAFF Capture Beads B7, 13X

Catalog# / Size	740549 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	TNFSF13B, BAFF, BLYS, CD257, tumor necrosis factor superfamily member 13b, DTL, TALL-1, TALL1, THANK, TNFSF20, ZTNF4, TNLG7A
Description	LEGENDplex™ Human B Cell Panel Capture Beads are intended for use with the following reagents: 740536 (LEGENDplex™ Human B Cell Panel Detection Antibodies) 740537 (LEGENDplex™ Human B Cell Panel Standard) 741150 (LEGENDplex™ Buffer Set A7)

Product Details

Reactivity	Human
Application	Multiplex Learn more about LEGENDplex™ at biolegend.com/legendplex Download the LEGENDplex™ software here .

Antigen Details

Biology Area	Cell Biology, Immunology, Signal Transduction, Stem Cells
Molecular Family	CD Molecules, Cytokines/Chemokines, Growth Factors, Soluble Receptors
Gene ID	10673

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587