

## LEGENDplex™ Human Resistin Capture Bead B7, 13X

<b>Catalog# / Size</b>	740208 / 100 tests (270 µL)
<b>Regulatory Status</b>	RUO
<b>Description</b>	This Human Resistin Capture Bead is suitable for use with the following Mix and Match panels:  <u>Human Adipokine Panel</u> 740211 (LEGENDplex™ Human Adipokine Detection Antibodies) 740210 (LEGENDplex™ Human Adipokine Standard) 740374 (LEGENDplex™ Buffer Set C)  <u>Human Metabolic Panel</u> 740214 (LEGENDplex™ Human Metabolic Panel 1 Detection Antibodies) 740213 (LEGENDplex™ Human Metabolic Panel 1 Standard) 740375 (LEGENDplex™ Buffer Set D)  For product manuals for these panels, please visit the Mix and Match tab on the <a href="#">LEGENDplex™</a> page.

### Product Details

---

<b>Reactivity</b>	Human
<b>Application</b>	<a href="#">Multiplex</a>  Learn more about LEGENDplex™ at <a href="http://biolegend.com/legendplex">biolegend.com/legendplex</a>  Download the <a href="#">LEGENDplex™ software here</a> .

### Antigen Details

---

<b>Molecular Family</b>	Enzymes and Regulators
<b>Gene ID</b>	<a href="#">56729</a>

### Related Protocols

---

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

