

LEGENDplex™ Mouse IgA Capture Bead B2, 13X

Catalog# / Size	740500 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	Ig Isotyping, IgG1, IgG2a, IgG2b, IgG3, IgA, IgM
Description	<p>LEGENDplex™ Mouse Immunoglobulin Isotyping Panel Capture Beads are intended for use with the following reagents:</p> <p>740494 (LEGENDplex™ Mouse Immunoglobulin Isotyping Panel Standard) 740495 (LEGENDplex™ Mouse Immunoglobulin Isotyping Panel Detection Antibodies) 740375 (LEGENDplex™ Buffer Set D)</p>

The LEGENDplex™ Mouse Immunoglobulin Isotyping Panel is a multiplex bead-based assay, using fluorescence–encoded beads suitable for use on various flow cytometers. This panel allows simultaneous quantification of 6 mouse proteins, including IgG1, IgG2a, IgG2b, IgG3, IgA, IgM. This assay panel provides high detection sensitivities and broad dynamic ranges. The panel has been validated by detecting expected changes in biological samples.

Product Details

Reactivity	Mouse
Application	Multiplex
	Learn more about LEGENDplex™ at biolegend.com/legendplex
	Download the LEGENDplex™ software here .

Antigen Details

Gene ID	NA
----------------	----

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587