

## PerCP/Cyanine5.5 anti-mouse CD19 Antibody

<b>Catalog# / Size</b>	152405 / 25 µg 152406 / 100 µg
<b>Clone</b>	1D3/CD19
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	B4, CVID3
<b>Isotype</b>	Rat IgG2a, κ
<b>Description</b>	CD19 is a 95 kD glycoprotein also known as B4. It is a member of the Ig superfamily, expressed on all pro-B to mature B cells (during development) and follicular dendritic cells. Plasma cells do not express CD19. CD19, in association with CD21 and CD81, forms a molecular complex integral to B cell activation.

### Product Details

---

<b>Verified Reactivity</b>	Mouse
<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Rat
<b>Immunogen</b>	Transfected k562 cells expressing the extracellular domain of murine CD19.
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
<b>Preparation</b>	The antibody was purified by affinity chromatography and conjugated with PerCP/Cyanine5.5 under optimal conditions.
<b>Concentration</b>	0.2 mg/ml
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. <b>Do not freeze.</b>
<b>Application</b>	<a href="#">FC - Quality tested</a>
<b>Recommended Usage</b>	Each lot of this antibody is quality control tested by <a href="#">immunofluorescent staining with flow cytometric analysis</a> . For flow cytometric staining, the suggested use of this reagent is ≤0.25 µg per million cells in 100 µl volume. It is recommended that the reagent be titrated for optimal performance for each application.  * PerCP/Cyanine5.5 has a maximum absorption of 482 nm and a maximum emission of 690 nm.
<b>Application Notes</b>	Additional reported applications (for the relevant formats) include: blocking <sup>1,2</sup> and immunoprecipitation <sup>1</sup> .
<b>Additional Product Notes</b>	BioLegend is in the process of converting the name PerCP/Cy5.5 to PerCP/Cyanine5.5. The dye molecule remains the same, so you should expect the same quality and performance from our PerCP/Cyanine5.5 products. Contact <a href="#">Technical Service</a> if you have any questions.
<b>Application References</b>	1. Krop I, <i>et al.</i> 1996. <i>Eur. J. Immunol.</i> 26:238. (FC, Blocking, IP) 2. Abraham P, <i>et al.</i> 2014. <i>Clin. Exp. Immunol.</i> 175:181. (Blocking)
<b>(PubMed link indicates BioLegend citation)</b>	
<b>Product Citations</b>	1. Fujita Y <i>et al.</i> 2018. <i>Cell reports.</i> 24(12):3296-3311 . <a href="#">PubMed</a> 2. McGinty JW, <i>et al.</i> 2020. <i>Immunity.</i> 52(3):528-541. <a href="#">PubMed</a> 3. Montel-Hagen A, <i>et al.</i> 2020. <i>Cell Rep.</i> 33:108320. <a href="#">PubMed</a> 4. Chen L, <i>et al.</i> 2020. <i>J Exp Med.</i> 217:00:00. <a href="#">PubMed</a> 5. Gogoi H, <i>et al.</i> 2020. <i>Front Immunol.</i> 1.620833333. <a href="#">PubMed</a> 6. Si Y, <i>et al.</i> 2020. <i>Sci Adv.</i> 6:eaba0995. <a href="#">PubMed</a>
<b>RRID</b>	AB_2629814 (BioLegend Cat. No. 152405)

## Antigen Details

<b>Structure</b>	Ig superfamily, associates with CD21 and CD81, 95 kD.
<b>Distribution</b>	Pro-B cells to mature B cells (during development), follicular dendritic cells.
<b>Function</b>	Modulates B cell activation and differentiation.
<b>Ligand/Receptor</b>	CD21, CD81, Leu-13.
<b>Cell Type</b>	B cells, Dendritic cells
<b>Biology Area</b>	Costimulatory Molecules, Immunology
<b>Molecular Family</b>	CD Molecules
<b>Antigen References</b>	<ol style="list-style-type: none"> <li>1. Fearon DT. 1993. <i>Curr. Opin. Immunol.</i> 5:341.</li> <li>2. Krop I, et al. 1996. <i>Eur. J. Immunol.</i> 26:238.</li> <li>3. Krop I, et al. 1996. <i>J. Immunol.</i> 157:48.</li> <li>4. Tedder TF, et al. 1994. <i>Immunol. Today</i> 15:437.</li> </ol>
<b>Gene ID</b>	<a href="#">12478</a>


## Related Protocols

[Cell Surface Flow Cytometry Staining Protocol](#)

## Other Formats

Purified anti-mouse CD19, FITC anti-mouse CD19, PerCP/Cyanine5.5 anti-mouse CD19, PE anti-mouse CD19, APC anti-mouse CD19, APC/Cyanine7 anti-mouse CD19, Alexa Fluor® 700 anti-mouse CD19

## Product Data



C57BL/6 mouse splenocytes were stained with CD19 (clone 1D3/CD19) PerCP/Cyanine5.5 (filled histogram) or rat IgG2a, κ PerCP/Cyanine5.5 isotype control (open histogram).

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587

