

TotalSeq™-D0092 Mouse IgG2b, κ Isotype Ctrl Antibody

Catalog# / Size	400383 / 10 µg
Clone	MPC-11
Regulatory Status	RUO
Isotype	Mouse IgG2b, κ
Barcode Sequence	ATATGTATCACGCGA
Description	The MPC-11 immunoglobulin has unknown specificity. This antibody was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse, rat and human tissues.

Product Details

Antibody Type	Monoclonal
Host Species	Mouse
Formulation	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 1 mM EDTA
Preparation	The antibody was purified by chromatography and conjugated with TotalSeq™-D oligomer under optimal conditions.
Concentration	0.5 mg/mL
Storage & Handling	The antibody solution should be stored undiluted between 2°C and 8°C. Do not freeze.
Application	PG - Quality tested
Recommended Usage	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis and the oligomer sequence is confirmed by sequencing. TotalSeq™-D antibodies are compatible with Mission Bio's Tapestry Single-Cell Sequencing Platform for simultaneous detection of DNA and Protein.

To maximize performance, it is strongly recommended that the reagent be titrated for each application, and that you centrifuge the antibody dilution before adding to the cells at 14,000xg at 2 - 8°C for 10 minutes. Carefully pipette out the liquid avoiding the bottom of the tube and add to the cell suspension. For Proteogenomics analysis, the suggested starting amount of this reagent for titration is ≤ 1.0 µg per million cells in 100 µL volume. Refer to the corresponding TotalSeq™ protocol for specific staining instructions.

Buyer is solely responsible for determining whether Buyer has all intellectual property rights that are necessary for Buyer's intended uses of the BioLegend TotalSeq™ products. For example, for any technology platform Buyer uses with TotalSeq™, it is Buyer's sole responsibility to determine whether it has all necessary third party intellectual property rights to use that platform and TotalSeq™ with that platform.

Application Notes	Additional reported applications (for the relevant formats) include: Intracellular Flow Cytometry (ICFC), Immunocytochemistry (ICC), Immunohistochemistry (IHC), Immunoprecipitation (IP), Western Blotting (WB), Functional Assay (FA).
--------------------------	--

Additional Product Notes	TotalSeq™-D reagents are designed to profile protein expression at single cell level. The Mission Bio Tapestry platform and sequencer (e.g. Illumina analyzers) are required. Please contact technical support for more information, or visit biolegend.com/totalseq/single-cell-dna
---------------------------------	--

The barcode flanking sequences are CGAGATGACTACGCTACTCATGG (PCR handle), and GAGCCGATCTAGTATCTCAGT*C*G (capture sequence). * indicates a phosphorothioated bond, to prevent nuclease degradation.

View more applications data for this product in our [Application Technical Notes](#).

Application References

(PubMed link indicates

1. Smed-Sørensen A, *et al.* 2008. *Blood* 111:5037. (FA) [PubMed](#)
2. Podolin PL, *et al.* 2008. *J. Immunol.* 180:7989. (FC) [PubMed](#)

BioLegend citation)

Antigen Details

Gene ID NA

Related Protocols

[TotalSeq™-D with Mission Bio Tapestri® Single-Cell DNA + Protein Protocol](#)

Other Formats

APC Mouse IgG2b, κ Isotype Ctrl, Biotin Mouse IgG2b, κ Isotype Ctrl, FITC Mouse IgG2b, κ Isotype Ctrl, PE Mouse IgG2b, κ Isotype Ctrl, PE/Cyanine5 Mouse IgG2b, κ Isotype Ctrl, Purified Mouse IgG2b, κ Isotype Ctrl, APC/Cyanine7 Mouse IgG2b, κ Isotype Ctrl, PE/Cyanine7 Mouse IgG2b, κ Isotype Ctrl, Alexa Fluor® 647 Mouse IgG2b, κ Isotype Ctrl, Alexa Fluor® 488 Mouse IgG2b, κ Isotype Ctrl, Pacific Blue™ Mouse IgG2b, κ Isotype Ctrl, Alexa Fluor® 700 Mouse IgG2b, κ Isotype Ctrl, PerCP Mouse IgG2b, κ Isotype Ctrl, PerCP/Cyanine5.5 Mouse IgG2b, κ Isotype Ctrl, Brilliant Violet 421™ Mouse IgG2b, κ Isotype Ctrl, Brilliant Violet 570™ Mouse IgG2b, κ Isotype Ctrl, Brilliant Violet 510™ Mouse IgG2b, κ Isotype Ctrl, Ultra-LEAF™ Purified Mouse IgG2b, κ Isotype Ctrl, Brilliant Violet 605™ Mouse IgG2b, κ Isotype Ctrl, Brilliant Violet 650™ Mouse IgG2b, κ Isotype Ctrl, Brilliant Violet 711™ Mouse IgG2b, κ Isotype Ctrl, Brilliant Violet 785™ Mouse IgG2b, κ Isotype Ctrl, PE/Dazzle™ 594 Mouse IgG2b, κ Isotype Ctrl, Alexa Fluor® 594 Mouse IgG2b, κ Isotype Ctrl, APC/Fire™ 750 Mouse IgG2b, κ Isotype Ctrl, GolnVivo™ Purified Mouse IgG2b, κ Isotype Ctrl, TotalSeq™-A0092 Mouse IgG2b, κ isotype Ctrl, Go-ChIP-Grade™ Purified Mouse IgG2b, κ isotype Ctrl, TotalSeq™-B0092 Mouse IgG2b, κ isotype Ctrl, TotalSeq™-C0092 Mouse IgG2b, κ isotype Ctrl, TotalSeq™-D0092 Mouse IgG2b, κ Isotype Ctrl, GMP APC Mouse IgG2b, κ Isotype Control

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587