

Recombinant Mouse CXCL12 (SDF-1 β) (carrier-free)

Catalog# / Size	589802 / 10 μ g 589804 / 25 μ g 589806 / 100 μ g
Regulatory Status	RUO
Other Names	Pre-B-cell growth-stimulating factor (PBSF), stromal cell derived factor 1 isoform beta precursor (SDF-1 β)
Description	CXCL12 is commonly known as SDF-1, with two variants, SDF-1 α (89 amino acids) and SDF-1 β (93 amino acids) derived from alternative splicing. They have identical amino acid sequences, differing only by the addition of four amino acids at the carboxy-terminus of the β form. CXCL12 β is constitutively expressed by stromal cells. CXCL12 α and β are ubiquitously expressed, and it has been suggested that their activity is regulated through differential proteolytic processing. Proteases present in serum process CXCL12 α initially at the carboxy-terminus (through carboxypeptidase N) and subsequently at the amino-terminus (through CD26); while CXCL12 β is only processed at the amino-terminus. These proteolytic processes affect the chemokine's ability to bind to heparin and cells as well as its ability to stimulate proliferation and chemotaxis. CXCL12 α is inactivated by other enzymes such as metalloproteases, serine proteases, and leukocyte elastase. CXCL12/CXCR4 pair plays a key role in maintaining proliferation and survival of bone marrow stem cell subpopulations under stress conditions. Also, CXCL12 is highly expressed in bone marrow stromal stem cells and regulates growth and survival. CXCL12 promotes survival by inducing anti-apoptotic proteins, decreasing apoptotic proteins, and increasing transcription of cell survival genes. Data suggests that CXCL12 β mediates cells' survival by enhancing autophagy in bone marrow mesenchymal stem cells.

Product Details

Source	Mouse CXCL12 (SDF-1 β) amino acids Lys 22-Met 93 (Accession# NP_038683.1) was expressed in <i>E. coli</i> .
Molecular Mass	The 72 amino acid recombinant protein has a predicted molecular mass of approximately 8.5 kD. The DTT-reduced and non-reduced protein migrate at approximately 11 and 14 kD by SDS-PAGE respectively. The N-terminal amino acid is Lysine.
Purity	>98%, as determined by Coomassie stained SDS-PAGE.
Formulation	0.22 μ m filtered protein solution is in PBS.
Endotoxin Level	Less than 0.01 ng per μ g cytokine as determined by the LAL method.
Concentration	10 and 25 μ g sizes are bottled at 200 μ g/mL. 100 μ g size and larger sizes are lot-specific and bottled at the concentration indicated on the vial. To obtain lot-specific concentration, please enter the lot number in our Concentration and Expiration Lookup or Certificate of Analysis online tools.
Storage & Handling	Unopened vial can be stored between 2°C and 8°C for up to 2 weeks, at -20°C for up to six months, or at -70°C or colder until the expiration date. For maximum results, quick spin vial prior to opening. The protein can be aliquoted and stored at -20°C or colder. Stock solutions can also be prepared at 50 - 100 μ g/mL in appropriate sterile buffer, carrier protein such as 0.2 - 1% BSA or HSA can be added when preparing the stock solution. Aliquots can be stored between 2°C and 8°C for up to one week and stored at -20°C or colder for up to 3 months. Avoid repeated freeze/thaw cycles.
Activity	Bioactivity was measured by its ability to chemoattract THP-1 cells in a dose dependent manner.
Application	Bioassay
Application Notes	BioLegend carrier-free recombinant proteins provided in liquid format are shipped on blue-ice. Our comparison testing data indicates that when handled and stored as recommended, the liquid format has equal or better stability and shelf-life compared to commercially available lyophilized proteins after reconstitution. Our liquid proteins are verified in-house to maintain activity after shipping on blue ice and are backed by our 100% satisfaction guarantee . If you have any concerns, contact us at tech@biolegend.com .

Additional Product Notes

View more applications data for this product in our [Scientific Poster Library](#).


Antigen Details

Structure	Chemokine
Distribution	Expressed by bone marrow stromal cells, endothelial cells, epithelial cells, pericytes, dendritic cells, fibroblasts
Function	CXCL12 is a constitutive chemokine involved in lung, brain, joint, and intestine inflammation. CXCL12 is chemotactic for T, B, and cancer cells. CXCL12/CXCR4 pair plays a role in mobilization of CD34 ⁺ progenitors to peripheral blood. CXCL12 is important in morphogenesis, angiogenesis, and cancer metastasis.
Interaction	Endothelial cells, epithelial cells, T and B cells, monocytes
Ligand/Receptor	CXCR4
Cell Type	B cells, Dendritic cells, Tregs
Biology Area	Cell Adhesion, Cell Biology, Immunology, Innate Immunity
Molecular Family	Adhesion Molecules, Cytokines/Chemokines
Antigen References	<ol style="list-style-type: none">1. Nagasawa T, <i>et al.</i> 1994. <i>Proc. Natl. Acad. Sci. USA</i> 91:2305.2. Muller A, <i>et al.</i> 2001. <i>Nature</i> 410:50.3. Valenzuela-Fernandez A, <i>et al.</i> 2002. <i>J. Biol. Chem.</i> 277:15677.4. De La Luz Sierra M, <i>et al.</i> 2004. <i>Blood</i> 103:2452.5. Kortesidis A, <i>et al.</i> 2005. <i>Blood</i> 105:3793.6. Davis DA, <i>et al.</i> 2005. <i>Blood</i> 105:4561.7. Zlotnik A. 2008. <i>J. Pathol.</i> 215:211.8. Herberg S, <i>et al.</i> 2013 <i>PLoS One</i> 8:e58207.

Gene ID

[6387](#)

Product Data



THP-1 cells chemoattracted by mouse CXCL12 (SDF-1β).

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587