

LEGENDplex™ Human Apolipoprotein (Apo) E4 Capture Bead B5, 13X

Catalog# / Size	740464 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	Apolipoprotein E4, ApoE4
Description	LEGENDplex™ Human Apolipoprotein (Apo) Panel Capture Beads are intended for use with the following reagents: 740456 (LEGENDplex™ Human Apolipoprotein (Apo) Panel Detection Antibodies) 740455 (LEGENDplex™ Human Apolipoprotein (Apo) Panel Standard) 740375 (LEGENDplex™ Buffer Set D)

The assembled Human Apolipoprotein (Apo) Mix and Match Panel is a multiplex bead-based assay panel, using fluorescence-encoded beads suitable for use on various flow cytometers. This allows simultaneous quantification of up to 11 human Apo proteins, including Apo AI, All, B100, CII, CIII, D, E, E4, H, J and M. This assay panel provides broader dynamic ranges than traditional ELISA methods. These targets have been validated for use with serum, plasma and cell culture supernatant samples.

Product Details

Reactivity	Human
Application	Multiplex Learn more about LEGENDplex™ at biolegend.com/legendplex Download the LEGENDplex™ software here .

Antigen Details

Biology Area	Immunology
Molecular Family	Enzymes and Regulators
Gene ID	348

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587