

LEGENDplex™ Human IL-18 Capture Bead B5, 13X

Catalog# / Size	740112 / 100 tests (270 µL)
Regulatory Status	RUO
Description	<p>This Human IL-18 Capture Bead is suitable for use with the following Mix and Match panels:</p> <p><u>Human Cytokine Panel 2 Panel</u> 740116 (LEGENDplex™ Human Cytokine Panel 2 Detection Antibodies) 740399 (LEGENDplex™ Human Cytokine Panel 2 Standard) 741156 (LEGENDplex™ Buffer Set Q)</p> <p><u>Human Neuroinflammation Panel 1</u> 740806 (LEGENDplex™ Human Neuroinflammation Panel 1 Detection Antibodies) 740807 (LEGENDplex™ Human Neuroinflammation Panel 1 Standard) 740368 (LEGENDplex™ Buffer Set A)</p> <p>For product manuals for these panels, please visit the Mix and Match tab on the LEGENDplex™ page.</p>

Product Details

Reactivity	Human
Application	Multiplex Learn more about LEGENDplex™ at biolegend.com/legendplex Download the LEGENDplex™ software here .

Antigen Details

Biology Area	Cell Biology, Immunology, Neuroinflammation, Neuroscience
Molecular Family	Cytokines/Chemokines
Gene ID	3606

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

