

Filter Plate for LEGENDplex™ Assay

Catalog# / Size	740377 / 1 plate 740378 / 5 plates
Regulatory Status	RUO
Other Names	Multiplex assay, Flow cytometry, Filter plates
Description	The LEGENDplex™ 96 well filter plates are available for use with LEGENDplex™ panels, mix & match systems, and buffer sets.

Product Details

Application	Multiplex Learn more about LEGENDplex™ at biolegend.com/legendplex Download the LEGENDplex™ software here .
--------------------	---

Antigen Details

Gene ID	NA
----------------	----

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587