

DAB Chromogen Concentrate

Catalog# / Size	926506 / 10 mL 926507 / 100 mL
Regulatory Status	RUO
Other Names	Chromogen Detection
Description	DAB is a commonly used chromogen for immunohistochemical staining. In the presence of a peroxidase enzyme, DAB will produce a brown precipitate that is insoluble in alcohol. The slides may be dehydrated through alcohols to xylene and mounted with any commercially available mounting media.

Product Details

Storage & Handling	Store between 2°C and 8°C.
Application	IHC
Recommended Usage	Each lot of this product is quality tested by immunohistochemistry. The optimal working dilution should be determined for each specific assay condition.

Antigen Details

Gene ID	NA
----------------	----

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587