

FITC anti-mouse CD154 (CD40L) Antibody

Catalog# / Size	157005 / 25 µg 157006 / 100 µg
Clone	SA047C3
Regulatory Status	RUO
Other Names	CD40L, gp39, TRAP, T-BAM, TNFSF5, CD40 Ligand
Isotype	Rat IgG2a, κ
Description	CD154 has been known a 39 kD TNF superfamily member also known as CD40 ligand, gp39, T-BAM, TRAP, and Ly-62. CD154 is an accessory molecule expressed predominantly on activated CD4 ⁺ lymphocytes that bind CD40. CD154 plays an important role in T-B cell costimulation and Ig class switching.

Product Details

Verified Reactivity	Mouse
Antibody Type	Monoclonal
Host Species	Rat
Immunogen	mCD40L-transfectants
Formulation	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
Preparation	The antibody was purified by affinity chromatography and conjugated with FITC under optimal conditions.
Concentration	0.5 mg/ml
Storage & Handling	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. Do not freeze.
Application	FC - Quality tested
Recommended Usage	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis . For flow cytometric staining, the suggested use of this reagent is ≤ 0.5 µg per million cells in 100 µl volume. It is recommended that the reagent be titrated for optimal performance for each application.
Excitation Laser	Blue Laser (488 nm)
RRID	AB_2814085 (BioLegend Cat. No. 157005) AB_2814086 (BioLegend Cat. No. 157006)

Antigen Details

Structure	A 39 kD TNF superfamily member
Distribution	Surface activated CD4 ⁺ T cells, B cells, macrophages, dendritic cells, mast cells, NK cells, polymorphs, monocytes, and activated platelets.
Function	Survival signaling to CD40 ⁺ B cells and Ig class switching.
Ligand/Receptor	CD40
Cell Type	B cells, Basophils, Dendritic cells, Fibroblasts, Macrophages, Mast cells, NK cells, Platelets, T cells
Biology Area	Costimulatory Molecules, Immunology

Molecular Family CD Molecules

Gene ID [959](#)


Related Protocols

[Cell Surface Flow Cytometry Staining Protocol](#)

Other Formats

Purified anti-mouse CD154 (CD40L), PE/Cyanine7 anti-mouse CD154 (CD40L), PE anti-mouse CD154 (CD40L), FITC anti-mouse CD154 (CD40L), APC anti-mouse CD154 (CD40L) Antibody, APC/Cyanine7 anti-mouse CD154 (CD40L), APC/Fire™ 750 anti-mouse CD154 (CD40L), PE/Dazzle™ 594 anti-mouse CD154 (CD40L), KIRAVIA Blue 520™ anti-mouse CD154 (CD40L)

Product Data



Enriched BALB/c mouse splenic T cells were stimulated with PMA plus ionomycin for 6 hours and then stained with mouse CD154 (CD40L) (clone SA047C3) FITC (filled histogram) or rat IgG2a FITC isotype control (open histogram).

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587