

LEGENDplex™ Human IL-1β Capture Bead B6, 13X

Catalog# / Size	740669 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	Interleukin 1 Beta, IL-1 Beta, Catabolin, IL1F2, Interleukin-1 Beta, IL1-BETA, IL-1
Description	This Human IL-1β Capture Bead is suitable for use in the following Mix and Match Panels: LEGENDplex™ Human Diabetes Panel 740916 (LEGENDplex™ Human Diabetes Panel Detection Antibodies) 740915 (LEGENDplex™ Human Diabetes Panel Standard) 740368 (LEGENDplex™ Buffer Set A)

Product Details

Reactivity	Human
Application	Multiplex Learn more about LEGENDplex™ at biolegend.com/legendplex Download the LEGENDplex™ software here .

Antigen Details

Biology Area	Cell Biology, Immunology, Innate Immunity, Neuroinflammation, Neuroscience
Molecular Family	Cytokines/Chemokines
Gene ID	3553

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587