

LEGENDplex™ SARS-CoV-2 Nucleocapsid Capture Bead B4, 13X

Catalog# / Size	741137 / 100 tests (270 µL)
Regulatory Status	RUO
Description	<p>The LEGENDplex™ SARS-CoV-2 Serological IgG Mix and Match Beads are individual capture beads for 3 key COVID-19 coronavirus derived proteins spike 1 (S1), receptor binding domain (RBD) and nucleocapsid (N). The beads are recommended for use with each other and the following Mix and Match assay reagents:</p> <p>LEGENDplex™ SARS-CoV-2 Serological IgG Panel Mix and Match: 741133 LEGENDplex™ SARS-CoV-2 Serological IgG Panel Detection Abs OR 741140 LEGENDplex™ SARS-CoV-2 Serological IgA Panel Detection Abs 741134 LEGENDplex™ SARS-CoV-2 Serological Panel Standard 740375 LEGENDplex™ Buffer Set D 740377 LEGENDplex™ Filter Plate OR 740379 LEGENDplex™ V-bottom Plate</p>

Product Details

Verified Reactivity	SARS-CoV-2
Application	Multiplex Learn more about LEGENDplex™ at biolegend.com/legendplex Download the LEGENDplex™ software here .

Antigen Details

Gene ID	43740575
----------------	--------------------------

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587

